

# CONDUCTIVE GEL



# Product Introduction

The conductive gel is made of hypoallergenic Aqueous electromedical coupling agent: Purified water, Carbopol, Glycerin, Sodium, hydroxide, Triclosan, etc., which optimizes electrostimulation sessions with comfort and efficiency. It is chloride-free; non-irritating; non-greasy; odorless, quickly absorbed. Safe and simple to use, and suitable for all types of skin, all types of electrodes.



## Usage Scenarios



# Product Features



## High Conductivity

It is designed to improve the conductivity by using a unique and highly conductive solution that facilitates better skin-to-device contact. Replacing layers of air with its aqueous solution, the gel removes any barriers that may reduce the machine's effectiveness and has been specially designed to transmit a variety of signals. In addition to increasing your device's efficiency, it will also increase its lifespan.



## Water-based gel

Contains 98% water, it is greasy; odorless so it, and can be quickly absorbed. It doesn't leave a residue and is easily removed with soap and water. It's also non-toxic, hypoallergenic, and safe for all skin types- even sensitive skin!



## Multiple application

It is suitable for a wide variety of electromedical uses.

- \*all Cavitation and Ultrasonic Slimming
- \*Ultrasound Devices
- \*TENS
- \*Radiofrequency facial & body machines
- \*EMS (electrical muscle stimulation), Ab belts & muscle stimulators
- \*Microcurrent devices (Nuface)
- \*At-home baby doppler
- \*Laser hair removal
- \*Electrical facial massagers
- \*Light therapy devices (blue and red light)



## Safe and easy to use

It is skin-friendly formula. It is easy to apply and requires no special equipment or technical know-how. Its simple instructions can be followed by anyone in just minutes.



## Increase the adhesion of electrodes

The user testing report shows, it can increase 30% adhesion of the electrode pad, which increase the using times of the electrodes!



## Customized Package






The containers can be customized and packed in tube, spray bottle, pump bottle, bags.



## Customized Formula

The containers can be customized, it can be packed in tube, spray bottle, pump bottle, bags.

## Product Specifications

Serial Number	Item No	Bottle type	Volume	Package
01	TG-65ML		65ml/2.3oz	200pcs/carton, Carton dimension: 43cm*36cm*28cm, Gross weight: 14.5kgs
02	TG-90ML		90ml/3.2oz	160pcs/carton, Carton dimension: 43cm*36cm*28cm, Gross weight: 16kgs
03	TG-100ML		100ml/3.5oz	160pcs/carton, Carton dimension: 48cm*40cm*33cm, Gross weight: 17.5kgs
04	TG-150ML		150ml/5.3oz	160pcs/carton, Carton dimension: 48cm*40cm*33cm, Gross weight: 17kgs
05	TG-250ML		250ml/8.9oz	50pcs/carton, Carton dimension: 40cm*30cm*30cm Gross Weight:17kgs

## FAQ

**Q** What is the shelf life of the conductive gel?

**A** The shelf life of the conductive gel is 2 year unopened. If it is opened, please used it up within 12 months.

**Q** Do your products contain latex?

**A** All products manufactured by TOP-RANK Healthcare are latex free. We do not use latex of any kind in the electrodes we manufacture or in the packaging.

**Q** Can I get free samples for testing?

**A** Yes, free samples is available, please contact us to get free sample!

### Office

Room2008, 20/F,Block A,China ResourcesBuild-  
ing,1366 Qianjiang Road, Hangzhou P.R.C(310020)

Tel: +86(575)82912999/82915999

Fax: +86(575)82670000

### Factory

MaShan Village, Dongguan Street, Shangyu  
City, Zhejiang Province, P.R.CHINA. 312352

Tel: +86(575)82563666

Fax: +86(575)82563750